

(CCMAR BIO-IMAGING



About us



Training, assistance and services in Bio-Imaging

TEAM



Teresa Matias Correia

Coordinator & Group Leader



Carina Mónico

Microscopy Service Specialist



Rute Castelo Félix

Team Member | Sample Preparation



Sandra Silva

Team Member | Clearing

CONTACT

ccmarimaging@ualg.pt

How can we help?



CONSULTATION

Meet us to discuss your needs and optimal bio-imaging approach.



TEST

(if necessary)
book a test with us on a specific microscope.

LARGE PROJECT

SMALL PROJECT



TRAINING

Mandatory training session to learn how to operate the microscope.



WE DO IT

Only need a few images?
Book the microscope and we do the acquisition for you. No training required.



ASSISTED USE

Get hands-on experience with our assistance. The training is finished when the researcher can operate the microscope safely and correctly.




INDEPENDENT USER

Book and use the microscope with regular user privileges. We are available for additional training or assistance.

PROCESSING & ANALYSIS

We can show you how to open and convert raw data to generic file formats, basic process & analysis procedures. We can also help you prepare a few images for a publication.



COLLABORATION

In you need, for example, custom-building hardware or software, novel (or uncommon/untested) protocols for sample preparation, image acquisition, processing and/or analysis, we can discuss a scientific collaboration.

Equipment



ZEISS LIGHTSHEET Z.1

3D Fluorescence Microscope



WORKSTATION

Image Processing and Analysis

LSFM FACILITY



ZEISS AXIO ZOOM.V16

Large-field Stereo Fluorescence Microscope



ZEISS AXIOSCOPE 5

Upright Fluorescence Microscope

BIOTERIUM FACILITY



ZEISS PRIMOVERT

Inverted Microscope with a Colour Camera



ZEISS AXIO VERT.A1


Inverted Fluorescence Microscope

CELL CULTURE FACILITY


Applications

ZEISS LIGHTSHEET Z.1

3D imaging of whole living or cleared samples.




Senegalese sole (flatfish)



Trout embryo

ZEISS AXIO ZOOM.V16


Observing, preparing, sorting and microsampling (dissecting, stimulating, manipulating) samples.



Zebrafish larvae blood vessels expressing fli1a:GFP

ZEISS AXIOSCOPE 5


Viewing samples with multiple fluorescent labels at high magnification (up to 100x).



BPAE cells stained to show nuclei (blue) and actin (green)

ZEISS PRIMOVERT


Cell culture routine work. Examine cell morphology and viability assays in brightfield and phase contrast.



Micrasterias radiata (algae), brightfield contrast

ZEISS AXIO VERT.A1

Investigate cell cultures using brightfield, phase contrast and multicolour fluorescence.



A549 cancer cells stained to show nuclei (blue), alpha-tubulin (green) and actin (red)

Zeiss Lightsheet Z.1

Building 8 | Lab 1.16



The Zeiss Lightsheet Z.1 is a lightsheet fluorescence microscope that enables fast volumetric imaging of large cleared samples (up to 1cm x 1cm x 2cm) or live transparent samples with low photobleaching and low phototoxicity.

LASER LINES

488 nm 50mW

561 nm 50mW

638 nm 75mW

OBJECTIVES | ILLUMINATION

5x/0.1

10x/0.2

OBJECTIVES | DETECTION

5x/0.16 EC Plan-NEOFLUAR

10x/0.5 W Plan-APO

20x/1.0 W Plan-APO DIC Corr UV-VIS-IR nd=1.33

20x/1.0 Clr Plan-APO Corr VIS-IR nd=1.38

FILTERS | CAMERA 1 & CAMERA 2

GFP | 505-545 nm & Cy3 | 575-615nm

GFP | 505-545 nm & mCherry | 585nm

YFP | 525-565 nm & Cy3 | 575-615nm

Cy3 | 575-615 nm & DRAQ5 | LP 660 nm

CAMERAS

2 PCO.edge sCMOS Cameras | 30 fps | pixel 6.5 μ m

CHAMBERS AND INCUBATION

Aqueous media (n=1.33)

Clearing media (n=1.38/20x, n=1.45/5x)

Incubation and temperature control options

ACQUISITION DESKTOP ZEN 2014


CPU 2x Intel XEON CPU E5-2623 @3.0GHz

RAM 64 GB

Storage 4x 2TB | 4TB RAID 10

GPU NIVIDIA Quadro K2000 4 GB

OS Windows 7




Zeiss Lightsheet Z.1 Workstation

Building 8 | Lab 1.16



Handling, processing and visualisation of lighsheet microscopy data is extremely computer intensive and requires the use of a high performance computer. The Zeiss Lightsheet Z.1 Win 7 workstation is specifically configured for storing and handling large experimental dataset obtained with the Zeiss Lightsheet Z.1 microscope.



WORKSTATION

CPU Intel XEON CPU E5-2620 @2.4GHz

RAM 192GB

Storage 6x HDD 8TB | RAID 5 config 36TB

GPU NVIDIA Quadro K2000 4 GB

OS Windows 7



SOFTWARE

Zeiss Zen 2014

Arivis vision4D

Fiji/ImageJ

WARNING

Save your data in D:\SWAP\your name and not in the C:\ drive.



Zeiss Axio Zoom.V16

Building 9 | Lab 0.08



The Zeiss Axio Zoom.V16 is a motorised fluorescence stereo zoom microscope that combines a large field of view, zoom and working distance. It is suitable for applications in which observation and careful handling of specimens exceeds the need for greater magnification. Applications include microinjections, dissections, movement or separation of small organisms, and more.

FEATURES

16x motorised zoom with high NA

Eyepieces 10x

Magnifications from 7x to 253x

Imaging modalities | brightfield, darkfield, oblique light, Epi-fluorescence

Acquisition of multi-channel images and **Z-stack**



ILLUMINATION

LED transmitted light for brightfield

Metal halide fluorescence lamp with 3 **fluorescence** excitation lines: UV (385) | Blue (475) | Green (555)



OBJECTIVES

1x/0.25 Plan-NEOFLUAR (WD 56mm)


2.3x/0.57 Plan-NEOFLUAR (WD 10.6mm)

FILTERS BLUE, GREEN, RED

49 DAPI | Ex 365 | BS 365 | EM BP 445/50

38HE GFP | Ex BP 470/40 | BS 495 | EM BP 525/50

63 mRFP | Ex BP 572/25 | BS 590 | EM BP 629/62



CAMERA

Axiocam 305 mono | 2.8 MP CMOS | 36 fps | pixel 4.54 μm

ACQUISITION DESKTOP ZEN 3.1 PRO

CPU Intel i7

RAM 8 GB

OS Windows 10

Zeiss Axioscope 5

Building 9 | Lab 0.08



The Zeiss Axioscope 5 is a manual upright microscope with transmitted white light LED and fluorescence LED illumination Colibri 7. The Zeiss Axioscope 5 has a wide selection of objectives, including a 100x oil lens for images with increased level of detail and resolution. It is a standalone microscope, which can acquire images without the need of a computer.

FEATURES

Imaging modalities | brightfield and epi-fluorescence
Eyepieces 10x



ILLUMINATION

Transmitted light LED 10W

Colibri 7 with 5 **fluorescence** excitation lines: UV (385) | Blue (475) | Green (555) | Yellow (590) | Red (630)

OBJECTIVES

10x/0.25 A-Plan

20x/0.45 A-Plan

40x/0.65 A-Plan

100x/1.25 A-Plan oil immersion



FILTERS BLUE, GREEN, RED

96 HE DAPI | Ex 390/40 | BS 420 | EM BP 450/40

38 HE GFP | Ex BP 470/40 | BS 495 | EM BP 525/50

45 HQ mCherry | Ex BP 560/40 | BS 585 | EM BP 630/75

CAMERA

Axiocam 202 mono | 2 MP CMOS | pixel 5.86 μ m



ACQUISITION DESKTOP ZEN 3.1 PRO

CPU Intel i7

RAM 8 GB

OS Windows 10

Zeiss Primovert

Building 7 | Lab 3.38



The Zeiss Primovert microscope is a manual inverted microscope equipped with the Axiocam 208 colour camera, designed for cell culture routine work.

FEATURES

Imaging modalities | brightfield and phase contrast
Eyepieces 10x



ILLUMINATION

Transmitted light HAL 30W

OBJECTIVES

4x/0.10 Plan-Achromat Ph0
10x/0.25 Plan-Achromat Ph1
20x/0.30 Plan-Achromat Ph1 LD
40x/0.50 Plan-Achromat Ph1 LD



PHASE RINGS

Ph0 and Ph1

CAMERA

Axiocam 105 colour | 5 MP CMOS | pixel 2.2 μm



BRING YOUR LAPTOP OR TABLET

Download Zen Lite available for free from the Zeiss Portal

Zeiss Axio Vert.A1

Building 7 | Lab 3.38



The Zeiss Axio Vert.A1 is a manual inverted microscope suitable for brightfield, phase contrast and epifluorescence to investigate cell cultures.

FEATURES

Imaging modalities | brightfield, phase contrast, and epifluorescence

Eyepieces 10x

Sample holders: suitable for multiwell plates, flasks, Petri dishes of 95 mm, microscope slides 76x26 mm, Ibidi chambers μ -SlideTM and Lab-TekTM II



ILLUMINATION

Transmitted light LED

Colibri 7 with 5 **fluorescence** excitation lines: UV (385) | Blue (475) | Green (555) | Yellow (590) | Red (630)

OBJECTIVES

5x/0.15 LD A-Plan

10x/0.25 LD A-Plan

20x/0.35 LD A-Plan

40x/0.55 LD A-Plan



FILTERS BLUE, GREEN, RED

96 HE DAPI | Ex 390/40 | BS 420 | EM BP 450/40

38 HE GFP | Ex BP 470/40 | BS 495 | EM BP 525/50

45 HQ mCherry | Ex BP 560/40 | BS 585 | EM BP 630/75

CAMERA

Axiocam 202 mono | 2 MP CMOS | pixel 5.86 μ m



ACQUISITION DESKTOP ZEN 3.1 PRO

CPU Intel i7

RAM 8 GB

OS Windows 10

How to access our equipment?

1



Email us at **ccmarimaging@ualg.pt** to schedule a consultation meeting to discuss your needs, experimental strategies and select the most adequate microscope.

Depending on your needs, we may acquire (a few) images for you or train you to become an independent user

2



Once you become an **independent user**, you will have access to the **online calendar** and will be able to **book** CCMAR's Bio-Imaging microscopes.

Note: "Assisted use" sessions need to take place within 1 week of the training session. If an independent user does not use a microscope for more than 6 months or misuses the microscope, we will contact the user to assess if a refresher session is required.

How to store and transfer data?


Save your data in E:\Users\your name (or D:\SWAP\your name).

Never save data in the C:\ drive.

Use **cloud** services to transfer data (e.g. Dropbox , Google drive, OneDrive).

Do not use external hard drives to transfer your data.

Note: Files older than 1 month will be removed from the local storage.



Usage Fees

**SPECIAL
OFFER**

TRY FOR FREE
until 31/05/2023

Instrument charge rates per hour from 1st of June 2023

Imaging System	CCMAR Internal Users	Academic visitors & non-profit organisations	Industry
Zeiss Lightsheet Z.1	10.00€	20.00€	As per quote
Zeiss Axio Zoom.V16	1.30€	2.60€	10.00€
Zeiss Axioscope 5	0.90€	1.80€	6.00€
Zeiss Axio Vert.A1	0.90€	1.80€	6.00€
Zeiss Primovert	0.70€	1.40€	4.00€

Usage Fees

CCMAR Bio-Imaging Services			
Services	CCMAR Internal Users	Academic visitors & non-profit organisations	Industry
Project Consultation and Advice	Free		As per quote
Training	Free	Free but instrument time is charged	
Testing Samples	Free	Free but instrument time is charged	
Sample Preparation	As per quote collaboration can be considered		
Assisted Imaging	As per quote collaboration can be considered		
Custom Image Analysis	As per quote collaboration can be considered		
Writing Support Imaging Methods (publications)	Free		
Writing Support Grant Applications (imaging sections)	Free		
Quotations for grant applications (imaging projects)	Free		
Advice concerning the purchase of microscopes and accessories	Free		
Microscope Maintenance and Quality Checks	As per quote		

Resources

STEP-BY-STEP INSTRUCTIONS

Microscope manuals coming soon!

ZEISS MANUALS

Microscope, sample preparation and image processing manuals

Zeiss Lightsheet Z.1 [sample preparation](#)

Zeiss Lightsheet Z.1 advances in [tissues clearing](#)

Zeiss Axio Zoom.V16 [info](#) and [manual](#)

Zeiss Axioscope 5 [info](#) and [manual](#)

ONLINE RESOURCES

Educational websites and videos

Zeiss webinars and how-to videos | [Zeiss Microscopy Insights Hub](#)

Zeiss Microscopy videos | [Zeiss Microscopy YouTube Channel](#)

Zeiss Lightsheet Z.1 Sample Preparation | [Zeiss How-to Video](#)

Zeiss Axioscope 5 Fluorescence Imaging | [Zeiss How-to Video](#)

Zeiss Axioscope 5 Transmitted Light Imaging | [Zeiss How-to Video](#)

Zeiss Education in Microscopy | [Carl Zeiss Microscopy Online Campus](#)

iBiology Microscopy Videos | [iBiology Microscopy Series](#)

Royal Microscopy Society Microscopy Videos | [RMS Technical Tea Breaks](#)

BIO-IMAGING SOFTWARE

Image visualization, processing and analysis software

Zeiss software downloads | [Zeiss Portal](#)

Zeiss Image Viewer | [Zeiss Zen Lite](#)

Open Source Image Processing Tool | [Fiji](#)

Fluorescence Spectra Viewer | [ThermoFisher](#)

Fluorescent Protein Database | [FPbase](#)

Grants, Publications & Acknowledgments

HELP WITH GRANT PROPOSALS

We can work with you on bio-imaging aspects of your grant application, help you write supporting statements and estimate microscopy access costs. We encourage you to discuss your planned work with us at the earliest opportunity.




MATERIALS & METHODS IN PUBLICATIONS

We can help you writing the Materials and Methods sections regarding microscopy and image processing.

Example: "Microscopy experiments were performed at the Bio-Imaging Facility at CCMAR."

"Images were acquired on a Zeiss [Axio Zoom.V16, Lightsheet Z.1 , etc] microscope equipped with a Zeiss [Axiocam 202 mono CMOS, etc] camera and controlled with the Zeiss Zen [3.1 Pro Blue Edition, etc] software, using the [1x and 0.25NA, etc] objective and the fluorescence filter sets [GFP and mCherry, etc].

"The [image analysis software] was used to generate [maximum intensity projections, 3D renderings, etc] of the acquired dataset.



ACKNOWLEDGMENTS IN PUBLICATIONS


When you use our equipment, please remember to acknowledge us and our funding sources in publications or other formal presentations. This is very important to guarantee future funding!



Acknowledgments: This study received Portuguese national funds from FCT - Foundation for Science and Technology through project UIDB/04326/2020, UIDP/04326/2020 and LA/P/0101/2020, and from the operational programmes CRESC Algarve 2020 and COMPETE 2020 through project EMBRC.PT ALG-01-0145-FEDER-022121. The authors would like to thank CCMAR Bio-Imaging for support & assistance in this work.

Collaborations: acknowledgments will be determined on a case by case basis.





CCMARIMAGING@UALG.PT

